

Student Name _____ Date _____

MASSACHUSETTS INSTITUTE OF TECHNOLOGY
Department of Electrical Engineering and Computer Science

6.161 Modern Optics Project Laboratory

Laboratory Exercise No. 4

Fall 2005

Holography

The Pre-Lab Exercises must be completed BEFORE entering the Lab. In your lab notebook record data, explain phenomena you observe, and answer the questions asked. Remember to answer all questions in your lab notebook in a neat and orderly fashion. No data are to be taken on these laboratory sheets. Tables provided herein are simply examples of how to record data into your laboratory notebooks. Expect the in-lab portion of this exercise to take about 3 hours.

PRE-LAB EXERCISES

PL4.1 – Get Prepared to Start the Laboratory Exercises

Read the **entire** laboratory handout, and be prepared to answer questions before, during and after the lab session. Determine all the equations and constants that may be needed in order to perform all the laboratory exercises. **Write** them all down in your laboratory notebook before entering the Lab. This will ensure that you take all necessary data while in the Lab in order to complete the lab write-up. This preparatory work will also count toward your Lab Exercise grade.

PL4.2 Preparation for Making Holograms

In Lab Exercise 4.2, your group will be making one type of hologram. You will have the option of making: (a) a Denisyuk single-beam white-light reflection hologram (simplest), (b) a two-beam or multiple beam reflection hologram, (c) a two-beam transmission hologram (must be viewed with a laser) or (d) a rainbow hologram (pseudo-color – but it takes quite a while to setup and develop). For those with little or no experience with optics or holography, this is a great opportunity to get your feet wet (pun intended).

For this section of the Pre-Lab it is your job to find out more about these types of holograms - everything about exposure of the holographic film, its development, and processing will be explained on lab day. However, before coming to the Lab you are advised to do some outside research on how to configure the holographic recording and readout setups for each of the four types of holograms listed above.

After you have decided which hologram you want to make, make a detailed sketch of the recording system you will use. – this includes placement of the object, lenses, laser, mirrors, beam splitters, holographic plate, and any other items you might need (e.g., slits, wave plates, irises, etc...).

Here are some useful tips to consider before coming to the MOL to make holograms:

6.161 - Laboratory Exercise No. 5, Fall 2005

- (1) You should bring some *interesting* objects from which we can make holograms. Such objects should not be dark in color, nor should they exceed a volume of 3"x3"x3".
- (2) Those students whose objects we use, will be able to take home a copy of their hologram. Otherwise, holograms made in the MOL will be archived for demonstration to future students. You will need your own laser to view the transmission holograms at home. Reflection holograms will be viewable in white light.
- (3) Make sure to wear clothes that can handle stains. Also, you will not be able to leave to go to the bathroom while the plates are developing, so be sure to go beforehand.
- (4) Be aware that it will be extremely dark at times in the laboratory, and that it can be fairly time-consuming depending on the complexity of the hologram you choose to make.

"By-the bye, what became of the baby?' said the Cat. 'I'd nearly forgotten to ask.'
'It turned into a pig,' Alice answered."

– Lewis Carroll, *Alice in Wonderland*

IN-LAB EXERCISES

Holography: An Example of Diffraction

4.1 Viewing of Transmission, Reflection, White-Light and Computer-Generated Holograms

A collimated laser beam and a white-light source have been set up for viewing of four types of holograms.

- (a) Illuminate the given laser-made **transmission holograms** with the divergent laser beam (remove the large collimating lens). For each hologram, find the direction of propagation of the virtual image, and view it directly (i.e., no lens or screen is to be used to view this image). Move your head around and tilt the hologram to get the best view of the image.

For the **hologram of the clock or the hologram of the lion**, draw a diagram of the setup that gives a virtual image with the best fidelity. In your diagram, show the location of the virtual image and any other beams that exit the hologram. Use this information to infer the geometry of the setup that was used to record the hologram (reverse engineering).

- (b) Illuminate the **reflection hologram of the coins** with the divergent laser beam and view the image through the hologram in reflection. Is the image you see real or virtual? Draw a diagram of the setup showing the location of the image and any other beams that exit the hologram. Use this information to infer the geometry of the setup used to record the hologram.
- (c) Illuminate the **white-light hologram of the owl** with white-light from the overhead lights in the Laboratory such that the light strikes the hologram and then reflects an image into your eye. Tip the hologram from side to side and up and down. Draw a diagram of your readout configuration showing the location of the image or images. Describe your observations from tipping the hologram in the beam. Re-illuminate the white light hologram with the divergent He-Ne laser beam. What differences do you observe? Explain these differences. Use your observations to infer the geometry of the setup used to record this hologram.
- (c) Put the collimating lens back into the system so as to generate a collimated laser beam. Illuminate the **first computer-generated transmission hologram** (written on a 4"x5" plate) with the collimated laser beam, and view the real image (it will spell a familiar acronym) on a distant screen. Draw a diagram of the setup showing the location of the real image and any other beams that exit the hologram. Use this information to infer the geometry of the setup used to record the hologram.
- (d) Place the **second computer-generated hologram** (written on 2"x2" film) into the collimated laser beam. Now relocate the screen close to second computer-generated hologram. Move the screen away from the hologram slowly. What do you observe on the screen? How do you explain what you see on the screen with respect to what you see on the actual hologram? Ask the TA or LA to show you how this hologram was made.

4.2 Making Your Hologram

In Prelab 4.2, you studied the recording setup for the following four types of holograms: (a) Denisyuk single-beam white-light reflection hologram, (b) two-beam or multiple beam reflection hologram, (c) two-beam transmission hologram and (d) rainbow hologram. Working with your group and with the assistance of the TA, now make one or more copies of the hologram type you chose.

- (a) In your lab notebook, describe and draw a diagram of the setup you used to make your hologram. Measure all distances between the optical and mechanical components in your setup, and show these dimensions in your diagram.
- (b) If possible, also specify: (a) the laser power density in the object and reference beams at the recording medium, and (b) the type, model number and manufacturer of the recording medium as appropriate.
- (c) Comment on the quality of the hologram that you made, and explain what could or should have been done to improve its quality.
- (d) Was this experience worthwhile?